



**2022**  
**GYNECOLOGIC**  
**BARIATRIC**  
**GENERAL**  
U.S. PRODUCT CATALOG

**LSI SOLUTIONS®**



**ADVANCE**  
SURGERY

---

**IMPROVE**  
OUTCOMES

---

**REDUCE**  
PAIN



Our approach at LSI SOLUTIONS® is that “Good Enough” is not good enough when helping a sick patient. We recognize that it is our privilege to truly enable a positive impact in a person’s life. We are dedicated to improving intra-operative and post-operative patient outcomes by advancing less invasive gynecologic surgery.

We believe surgery is now poised for an amazing Renaissance. This new era must be kinder and gentler. Surgery must focus on offering patients highly effective durable clinical outcomes with minimal pain and rapid recovery. We are proud to encourage, equip, and empower surgeons striving for optimized patient care. LSI will help lead this important advance.

Our goal is to help more patients in more places through our growing portfolio of exciting products. LSI is dedicated to supporting wonderful surgeons and their teams across the world with innovative technology, techniques and training. We are grateful to provide meaningful value to society, and most importantly, to patients.

A handwritten signature in black ink that reads "Jude" with a stylized flourish underneath.

Jude S. Sauer, MD  
Founder, President and CEO

# Contents

## Suturing

4

### **RD180° DEVICE 5MM**

Automated Laparoscopic Suturing Device - 5mm Shaft Diameter - 33cm Length

6

### **RD180° DEVICE 10MM**

Automated Laparoscopic Suturing Device - 10mm Shaft Diameter - 33cm Length

6

### **RD180° SP DEVICE**

Automated SINGLE PORT Suturing Device - 6.35mm Shaft Diameter - 43cm Length

6

### **RD° THE RUNNING DEVICE®**

Automated Suturing Device - 5mm Shaft Diameter - 33cm Length

6

### **RD° BARIATRIC LENGTH DEVICE**

Automated Suturing Device - 5mm Shaft Diameter - 43cm Length

7

---

### **RD° QUICK LOAD® SUTURE 2-0 STRONGSORB®**

2-0 PGA, Absorbable, 53", Violet, Braided

7

### **RD° QUICK LOAD® SUTURE 2-0 MONOGLIDE®**

2-0 PDO, Absorbable, 53", Purple, Monofilament

7

### **RD° QUICK LOAD® SUTURE 0 MONOGLIDE®**

0 PDO, Absorbable, 53", Purple, Monofilament

7

### **RD° QUICK LOAD® SUTURE 2-0 POLYESTER**

2-0 Polyester, Non-Absorbable, 53", Green, Braided

8

### **RD° QUICK LOAD® SUTURE 0 POLYESTER**

0 Polyester, Non-Absorbable, 53", Green, Braided

8

### **RD° QUICK LOAD® SUTURE 2-0 POLYPROPYLENE**

2-0 Polypropylene, Non-Absorbable, 53", Blue, Monofilament

8

## Knotting

10

### **TK° TI-KNOT® DEVICE**

Titanium Suture Fastening Device - 5mm Shaft Diameter - 31cm Length

12

### **TK° QUICK LOAD® UNIT**

Titanium Fastener in Customized Quick Load Unit



12

<i>multiSTATION™ SYSTEM triDOC™ Rail Clamp</i>	16
<i>multiSTATION™ SYSTEM duoDOC™ L Rail Clamp</i>	16
<i>multiSTATION™ SYSTEM duoDOC™ R Rail Clamp</i>	16
<i>multiSTATION™ SYSTEM monoDOC™ Rail Clamp</i>	16
<i>multiSTATION™ SYSTEM M45™ CLICK RELEASE™ miniARM® HOLDER</i>	17
<i>multiSTATION™ SYSTEM M55™ RIGID miniARM® HOLDER</i>	17
<i>multiSTATION™ SYSTEM SubX™ STERNAL ASCENDER™ SYSTEM</i>	17
<i>multiSTATION™ SYSTEM LIMA RetroSterno™ Paddle</i>	17
<i>multiSTATION™ SYSTEM RIMA RetroSterno™ Paddle</i>	18
<i>multiSTATION™ SYSTEM RetroSterno™ Paddle Kit</i>	18
<i>multiSTATION™ SYSTEM miniARM® LOCKING COUPLER</i>	18
<i>multiSTATION™ SYSTEM miniARM® CANNULA KIT</i>	18
<i>multiSTATION™ SYSTEM miniARM® CLAMPING ADAPTER</i>	19
<i>multiSTATION™ SYSTEM miniARM® ACCESS ADAPTER</i>	19
<i>multiSTATION™ SYSTEM miniARM® STERILIZATION TRAY</i>	19
<hr/>	
<b>FS® SOUND 6cm</b>	
Reusable 6cm Sound with Rotational Anchor	20
<b>FS® SOUND 8cm</b>	
Reusable 8cm Sound with Rotational Anchor	20
<b>FS® SOUND 10cm</b>	
Reusable 10cm Sound with Rotational Anchor	20
<b>FS® SOUND 12cm</b>	
Reusable 12cm Sound with Rotational Anchor	20
<b>FS® SOUND SET</b>	
Reusable Sounds with Rotational Anchors - Set	21
<b>FS-30 FORNISEE® DEVICE</b>	
Illuminated Uterine Manipulator - 30mm Cervical Cup	21
<b>FS-35 FORNISEE® DEVICE</b>	
Illuminated Uterine Manipulator - 35mm Cervical Cup	21
<b>FS-40 FORNISEE® DEVICE</b>	
Illuminated Uterine Manipulator - 40mm Cervical Cup	21
<b>FS® DEVICE COMBO BOX</b>	
Illuminated Uterine Manipulators - Set	22
<b>FS® LIGHT GUIDE</b>	
Fiber Optic Illumination Cable	22
<hr/>	
<b>PNEUMOSTOP® DEVICE</b>	
Vaginal Canal Occluder	22
<b>PNEUMOSTOP® PLUS DEVICE</b>	
Vaginal Canal Occluder	22

**All products in this catalog are supplied sterile unless otherwise specified.**

# SUTURING

We strive to move the needle forward. Our automated suturing technology is designed to enable minimally invasive procedures, reduce trauma, save time and make challenging anatomy easier for you and your team to sew.



# Suturing

## **RD180® DEVICE 5MM**

Automated Laparoscopic Suturing Device - 5mm Shaft Diameter – 33cm Length



Reorder	Description	Quantity
021100	<b>RD180® DEVICE 5MM</b> 1 Device per Box	Box of 6 Devices

## **RD180® DEVICE 10MM**

Automated Laparoscopic Suturing Device - 10mm Shaft Diameter – 33cm Length



Reorder	Description	Quantity
021330	<b>RD180® DEVICE 10MM</b> 1 Device per Box	Box of 6 Devices

## **RD180® SP DEVICE**

Automated SINGLE PORT Suturing Device - 6.35mm Angled Shaft Diameter – 43cm Length



Reorder	Description	Quantity
021450	<b>RD180® SP DEVICE</b> 1 Device per Box	Box of 6 Devices

## **RD® THE RUNNING DEVICE®**

Automated Suturing Device - 5mm Shaft Diameter – 33cm Length



Reorder	Description	Quantity
021000	<b>RD® THE RUNNING DEVICE®</b> 1 Device per Box	Box of 6 Devices



**RD® THE RUNNING DEVICE® BARIATRIC**

Automated Suturing Device - 5mm Shaft Diameter – 43cm Bariatric Length



Reorder	Description	Quantity
021021	<b>RD® THE RUNNING DEVICE® BARIATRIC</b> 1 Device per Box	Box of 6 Devices

**RD® QUICK LOAD® SUTURE 2-0 STRONGSORB®**

2-0 PGA, Absorbable, 53", Violet, Braided



Reorder	Description	Quantity
021025	<b>RD® QUICK LOAD® SUTURE 2-0 STRONGSORB®</b> 2-0 PGA, Absorbable, 53", Violet, Braided 1 Suture per Pouch	Box of 12 Sutures

**RD® QUICK LOAD® SUTURE 2-0 MONOGLIDE®**

2-0 PDO, Absorbable, 53", Purple, Monofilament



Reorder	Description	Quantity
021030	<b>RD® QUICK LOAD® SUTURE 2-0 MONOGLIDE®</b> 2-0 PDO, Absorbable, 53", Purple, Monofilament 1 Suture per Pouch	Box of 12 Sutures

**RD® QUICK LOAD® SUTURE 0 MONOGLIDE®**

0 PDO, Absorbable, 53", Purple, Monofilament



Reorder	Description	Quantity
021506	<b>RD® QUICK LOAD® SUTURE 0 MONOGLIDE®</b> 0 PDO, Absorbable, 53", Purple, Monofilament 1 Suture per Pouch	Box of 12 Sutures

# Suturing

---

## **RD® QUICK LOAD® SUTURE 2-0 POLYESTER**

2-0 Polyester, Non-Absorbable, 53", Green, Braided



Reorder	Description	Quantity
020979	<b>RD® QUICK LOAD® SUTURE 2-0 POLYESTER</b> 2-0 Polyester, Non-Absorbable, 53", Green, Braided 1 Suture per Pouch	Box of 12 Sutures

---

## **RD® QUICK LOAD® SUTURE 0 POLYESTER**

0 Polyester, Non-Absorbable, 53", Green, Braided



Reorder	Description	Quantity
021010	<b>RD® QUICK LOAD® SUTURE 0 POLYESTER</b> 0 Polyester, Non-Absorbable, 53", Green, Braided 1 Suture per Pouch	Box of 12 Sutures


---

## **RD® QUICK LOAD® SUTURE 2-0 POLYPROPYLENE**

2-0 Polypropylene, Non-Absorbable, 53", Blue, Monofilament



Reorder	Description	Quantity
020989	<b>RD® QUICK LOAD® SUTURE 2-0 POLYPROPYLENE</b> 2-0 Polypropylene, Non-Absorbable, 53", Blue, Monofilament 1 Suture per Pouch	Box of 12 Sutures



# KNOTTING

Stronger, faster, and more reliable than hand-tied knots. Our titanium suture fastening technology is advancing surgery, improving outcomes, and reducing pain.



# Knotting

## **TK® TI-KNOT® DEVICE**

Titanium Suture Fastening Device  
5mm Shaft Diameter – 31cm Length



Reorder	Description	Quantity
030404	<b>TK® TI-KNOT® DEVICE</b> 1 Device per Box	Box of 6 Devices

## **TK® QUICK LOAD® UNIT SINGLES**

Titanium Fastener in Customized Quick Load Unit





Reorder	Description	Quantity
030510	<b>TK® QUICK LOAD® UNIT SINGLES</b> 1 Fastener per Pouch	Box of 12 Singles



# ACCESS & EXPOSURE

From Instrument holders to uterine manipulation, our technology is designed to assist, improve ergonomics, reduce trauma, and make surgery easier for you and your team.





# Access & Exposure

## multiSTATION™ SYSTEM

### triDOC™ Rail Clamp



Reorder	Description	Quantity
081000	multiSTATION™ SYSTEM triDOC™ Rail Clamp 1 Rail Clamp per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### duoDOC™ L Rail Clamp



Reorder	Description	Quantity
081023	multiSTATION™ SYSTEM duoDOC™ L Rail Clamp 1 Rail Clamp per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### duoDOC™ R Rail Clamp



Reorder	Description	Quantity
081024	multiSTATION™ SYSTEM duoDOC™ R Rail Clamp 1 Rail Clamp per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### monoDOC™ Rail Clamp



Reorder	Description	Quantity
081022	multiSTATION™ SYSTEM monoDOC™ Rail Clamp 1 Rail Clamp per Box Supplied Non-Sterile	Box of 1 Device

# Access & Exposure

## multiSTATION™ SYSTEM

### M45™ CLICK RELEASE™ miniARM® HOLDER



Reorder	Description	Quantity
081010	multiSTATION™ SYSTEM M45™ CLICK RELEASE™ miniARM® HOLDER 1 Holder per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### M55™ RIGID miniARM® HOLDER



Reorder	Description	Quantity
080950	multiSTATION™ SYSTEM M55™ RIGID miniARM® HOLDER 1 Holder per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### SubX™ STERNAL ASCENDER™ SYSTEM



Reorder	Description	Quantity
080900	multiSTATION™ SYSTEM SubX™ STERNAL ASCENDER™ SYSTEM 1 Sternal Ascender per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### LIMA RetroSterno™ Paddle



Reorder	Description	Quantity
081214	multiSTATION™ SYSTEM LIMA RetroSterno™ Paddle 1 Paddle per Box Supplied Non-Sterile	Box of 1 Device

# Access & Exposure

## multiSTATION™ SYSTEM

### RIMA RetroSterno™ Paddle



Reorder	Description	Quantity
081241	multiSTATION™ SYSTEM RIMA RetroSterno™ Paddle 1 Paddle per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### RetroSterno™ Paddle Kit



Reorder	Description	Quantity
081202	multiSTATION™ SYSTEM RetroSterno™ Paddle Kit 1 LIMA and 1 RIMA Paddle per Box Supplied Non-Sterile	Box of 1 Kit

## multiSTATION™ SYSTEM

### miniARM® LOCKING COUPLER



Reorder	Description	Quantity
080555	multiSTATION™ SYSTEM miniARM® LOCKING COUPLER For Endoeyeflex® 5mm 1 Locking Coupler per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

### miniARM® CANNULA ADAPTER KIT



Reorder	Description	Quantity
080415	multiSTATION™ SYSTEM miniARM® CANNULA ADAPTER KIT 1 Cannula Adapter, 1 Cannula 5mm x 17cm, and 1 Obturator per Box Supplied Non-Sterile	Box of 1 Kit

# Access & Exposure

## multiSTATION™ SYSTEM

miniARM® CLAMPING ADAPTER - Shafted Instrument Holder



Reorder	Description	Quantity
080519	multiSTATION™ SYSTEM miniARM® CLAMPING ADAPTER 1 Adapter per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

miniARM® ACCESS ADAPTER - 3D® MAP Holder



Reorder	Description	Quantity
080670	multiSTATION™ SYSTEM miniARM® ACCESS ADAPTER 1 Adapter per Box Supplied Non-Sterile	Box of 1 Device

## multiSTATION™ SYSTEM

miniARM® STERILIZATION TRAY®



Reorder	Description	Quantity
100034	multiSTATION™ SYSTEM miniARM® SYSTEM STERILIZATION TRAY Will accept all multiSTATION™ SYSTEM and miniARM® products 1 Tray per Box Supplied Non-Sterile	Box of 1 Device

Manufactured by: Summit Medical, St. Paul, MN 55121

# Access & Exposure

## **FS® SOUND 6CM**

Reusable Sound with Rotational Anchor – 6cm Angled Shaft



Reorder	Description	Quantity
110056	<b>FS® SOUND 6CM</b> 1 Device per Pouch Supplied Non-Sterile	1 Reusable Device

## **FS® SOUND 8CM**

Reusable Sound with Rotational Anchor – 8cm Angled Shaft



Reorder	Description	Quantity
110058	<b>FS® SOUND 8CM</b> 1 Device per Pouch Supplied Non-Sterile	1 Reusable Device

## **FS® SOUND 10CM**

Reusable Sound with Rotational Anchor – 10cm Angled Shaft



Reorder	Description	Quantity
110105	<b>FS® SOUND 10CM</b> 1 Device per Pouch Supplied Non-Sterile	1 Reusable Device

## **FS® SOUND 12CM**

Reusable Sound with Rotational Anchor – 12cm Angled Shaft



Reorder	Description	Quantity
110106	<b>FS® SOUND 12CM</b> 1 Device per Pouch Supplied Non-Sterile	1 Reusable Device

# Access & Exposure

## FS® SOUND SET

Reusable Sounds with Rotational Anchors – Set



Reorder	Description	Quantity
110112	<b>FS SOUND SET</b> Set of Four Reusable Sounds 6cm, 8cm, 10cm, 12cm Supplied Non-Sterile	4 Reusable Sounds

## FS-30 FORNISEE® DEVICE

Illuminated Uterine Manipulator – 30mm Cervical Cup



Reorder	Description	Quantity
110030	<b>FS-30 FORNISEE® DEVICE</b> 1 Device per Box	Box of 6 Devices

## FS-35 FORNISEE® DEVICE

Illuminated Uterine Manipulator – 35mm Cervical Cup



Reorder	Description	Quantity
110035	<b>FS-35 FORNISEE® DEVICE</b> 1 Device per Box	Box of 6 Devices

## FS-40 FORNISEE® DEVICE

Illuminated Uterine Manipulator – 40mm Cervical Cup



Reorder	Description	Quantity
110040	<b>FS-40 FORNISEE® DEVICE</b> 1 Device per Box	Box of 6 Devices

# Access & Exposure

## **FS® DEVICE COMBO BOX**

Illuminated Uterine Manipulators – Set



Reorder	Description	Quantity
<b>110029</b>	<b>FS® DEVICE COMBO BOX</b> Combo Box of 6 Single Patient Use Devices (2) FS-30 Devices (2) FS-35 Devices (2) FS-40 Devices	Combo Box of 6 Devices (2 of each size)

## **FS® LIGHT GUIDE**

Fiber Optic Illumination Cable



Reorder	Description	Quantity
<b>100025</b>	<b>FS® LIGHT GUIDE</b> 1 Reusable Device Supplied Non-Sterile	Box of 1 Device

Manufactured by: Fiberoptics Technology, Inc., Pomfret, CT 06258

## **PNEUMOSTOP® DEVICE**

Vaginal Canal Occluder



Reorder	Description	Quantity
<b>110200</b>	<b>PNEUMOSTOP® DEVICE</b> 3 Devices per Box	6 Boxes of 3 Devices

## **PNEUMOSTOP® PLUS DEVICE**

Vaginal Canal Occluder



Reorder	Description	Quantity
<b>110230</b>	<b>PNEUMOSTOP® PLUS DEVICE</b> 3 Devices per Box	6 Boxes of 3 Devices





**Photo Credits**

**Front Cover and Page 11:**

Tahir Chauhdry, MD, Doctor of Obstetrics, Olean General Hospital, Olean, NY

**Page 5:**

Steven McCarus, MD, Chief of Gynecological Surgery at AdventHealth Celebration





# INSPIRED BY **CARING**

Our customer is ultimately  
a patient. Our technology  
challenges human illness.  
Our mission is to lead the  
world in surgical innovation  
through product  
development and sales.



7796 VICTOR-MENDON RD | VICTOR, NY 14564 USA | CUSTOMER SERVICE: 866.575.3493 | TECH SUPPORT: 866.428.9092 | FAX: 585.742.8086 | [WWW.LSISOLUTIONS.COM](http://WWW.LSISOLUTIONS.COM)

The LSI logo, LSI SOLUTIONS, Click Release, duoDoc, Fornisee, FS, Heart logo, Inspired by Caring, M45, M55, miniARM, monoDoc, Monoglide, multiStation, Pneumostop, Quick Load, RD, RD180, RetroSterno, Sternal Ascender, Strongsorb, SubX, The Running Device, Ti-Knot, TK, and triDoc, are trademarks and registered trademarks of LSI SOLUTIONS, Inc. Copyright © 2020, LSI SOLUTIONS®. All Rights Reserved.

PN032455B 1.27.22