












RD180® PRODUCT ORDERING		SUPPLIED: STERILE INDIVIDUALLY PACKAGED	
	REORDER	PRODUCT	DESCRIPTION
 x 6	REF 021100	5mm RD180® DEVICE	Box of 6 Devices
 x 6	REF 021330	10mm RD180® DEVICE	Box of 6 Devices
 x 6	REF 021450	SP™ RD180® DEVICE	Box of 6 Devices
ABSORBABLE	 x 12	STRONGSORB® RD® QUICK LOAD® SUTURE 2-0 PGA, Braided, Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
	 x 12	MONOGLIDE® RD® QUICK LOAD® SUTURE 2-0 PDO, Monofilament, Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
	 x 12	MONOGLIDE® RD® QUICK LOAD® SUTURE 0 PDO, Monofilament, Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
NON-ABSORBABLE	 x 12	POLYESTER RD® QUICK LOAD® SUTURE 2-0 POLYESTER, Braided, Non-Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
	 x 12	POLYESTER RD® QUICK LOAD® SUTURE 0 POLYESTER, Braided, Non-Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)
	 x 12	POLYPROPYLENE RD® QUICK LOAD® SUTURE 2-0 POLYPROPYLENE, Monofilament, Non-Absorbable, 53"	Box of 12 Sutures (1 Suture per Pouch)

PN022880B 08.03.2022



CUFF CLOSURE | 5mm RD180® DEVICE



CUFF CLOSURE | 10mm RD180® DEVICE





RD180® DEVICE

Automated Laparoscopic Suturing Devices

A line of automated suturing devices specifically designed for laparoscopic cuff closure.



5mm RD180® DEVICE suturing the vaginal cuff.



NO EXPOSED SHARPS

The needle is always retracted and protected in distal shaft of device until lever is squeezed.

5mm, 10mm, OR SINGLE PORT SHAFT



Choose between three configurations for remote suture placement during minimally invasive procedures.

PINK LEVER

1 SQUEEZE® technology enables easy suture placement and suture rearm in a single squeeze of the lever.

ERGONOMIC HANDLE


Right or left hand comfort.



5mm | Allows for smaller ports and ideal for those patients interested in a more cosmetic hysterectomy



10mm | Larger device tip allows for extra large tissue bites incorporating mucosa and peritoneum in one bite



SP™ | Longer, angled shaft designed for ergonomics and enhanced visibility during single port hysterectomy

