


LONGER
ROTATING SHAFT

MORE REACH
EASIER CONTROL




LSOLUTIONS®

7796 Victor-Mendon Road
Victor, New York 14564 U.S.A.
Customer Service: +1.866.575.3493
Technical Support: +1.866.428.9092
Fax: +1.585.742.8086
www.lsisolutions.com

NEW
COR-KNOT
MINI® DEVICE


COR-KNOT MINI[®] DEVICE

SPINI MINI SHAFT | 17cm Long, Rotating Shaft




LS SOLUTIONS COR-KNOT MINI

SLIDE




Suture Slot Faces Center

SQUEEZE & HOLD



Compression & Tension

TUG & RELEASE



Squeeze Lever Only Once

REMOVE DEVICE








Inspect Fastener



COR-KNOT MINI[®] PRODUCT ORDERING

SUPPLIED: STERILE

	REORDER	PRODUCT	DESCRIPTION
 x 6	REF 031400	COR-KNOT MINI [®] DEVICE KIT	Box of 6 Kits (2 Devices per Kit)
 x 12	REF 030950	COR-KNOT [®] QUICK LOAD [®] SINGLES	Box of 12 Singles (1 Fastener per Pouch)
 x 12	REF 030902	COR-KNOT [®] QUICK LOAD [®] 6-POUCH	Box of 12 Pouches (6 Fasteners per Pouch)
 x 6  x 6	REF 031450	COR-KNOT MINI [®] COMBO KIT	Box of 6 Kits (2 Devices & 12 Fasteners per Kit)

PN031441A 06.15.16

LS SOLUTIONS[®]